

A dynamic background image featuring a splash of clear water against a teal gradient. The water is captured in mid-air, with numerous bubbles and droplets visible, creating a sense of movement and freshness. The teal gradient is darker at the bottom and lighter at the top, blending into the white background of the top left corner.

# *Webinar* Photometry & Automated Chemistry Analyzers

XAVIER TAN -&- BENJAMIN CHIANG

# WTW Webinar Series



**June 25**  
How to monitor BOD  
with OxiTop  
(Part #1)



**July 7**  
Online Wastewater  
Monitoring of COD/BOD  
With Spectral Sensors



**July 9**  
How to monitor BOD  
with OxiTop  
(Part #2)



**July 14**  
Photometry basics  
& Automated Chemistry  
Analyzers

# GoToWebinar

## Poll Question

QUICKPOLL

### What does Xylem mean?

Please select one:

- The transport tissue in plants which conducts water
- Chewing gum
- A global water technology provider

**Submit**

## Control Panel

The screenshot shows the GoToWebinar control panel. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below this is the 'Audio' section, which includes a 'Sound Check' indicator, radio buttons for 'Computer audio' (selected) and 'No audio', a 'MUTED' status, and dropdown menus for 'Microphone Array (Realtek(R) Audio)' and 'Speakers/Headphones (Realtek(R) Au...'. A volume slider is also present. Below the audio settings is the 'Questions' section, which is highlighted with a red bracket and the label 'Questions'. This section contains a 'Send' button and a text input field. At the bottom of the control panel, it displays 'test webinar', 'Webinar ID: 790-366-091', a recording status 'This session is being recorded.', and the GoToWebinar logo.

# WTW Photometry



# Xavier Tan

- Lab and Process Product Manager for water and wastewater, South East Asia
- 8 years laboratory and online process instrumentation
- 1 year with Xylem
- Product manager for WTW, STM & MJK instrumentation





**Poll Question #1**

Which part of the industry are you from?


# Overview

- Introduction to Photometry
  - Light & Optics
  - Measurements
  
- Xylem Analytical Lab & Field Photometric Solutions
  - WTW Lab & Field Photometers
  - OI Automated Chemistry Analyzer
  
- Q & A

# Introduction to Photometry

**Phos** (Greek) for light => Photometry is a measurement method to analyse (aqueous) solutions by means of a light source.

Light (physical) is a spectra of electromagnetic waves, divided into different ranges: Visible light (white light) ranges from approx. 380 – 780 nm



## **Photometric / Colorimetric Analysis:**

Determination of substances by their specific colour reaction / absorbance to their chemical properties at a specific wavelength.



# Introduction – Light and Optics


## Specific Wavelengths are obtained by

### Different light sources

- LEDs ( $\lambda_x$ ) = lowest power consumption, lower intensity
- Tungsten (White light halogene lamp for VIS range)
- Xenon (UV-VIS) => Flash lamp with „lifetime“ durability
- Deuterium (UV) => special lamp, expensive


### and different optical techniques

- Monochromator
- Polychromator
- Filter
- others




# Optics: Filter and LED Photometer

Filter photometer with reference beam: photoLab® S6/S12 Series




LED<sub>λ</sub> + optical Filter – single beam: photoFlex® Series



# Optics - Spectrophotometer

## Monochromator of photoLab® 7000 Series



### Optical Elements of photoLab:

**Lamp:** Halogenlamp (7100), Xenon-Flash lamp (7600)

**Mono-chromator:** This optical element is a optical band pass to uncouple light for required measuring wavelength

**Shutter:** Fades in and out the light beam to suppress ambient light (some photometer models only)

**Lenses:** Geometric beam formation for cuvette passing light beam

**Optical Filter** Reduction of straylight by additional supression of light of „wrong“ wavelength (advantage here: close to diode)

**Diode** Opto-electronic detector

# Measurement modes

**What type of measurement is performed in (UV-)VIS?**

**=> 3 Measurement modes of Photometric Analysis**


1.) **Transmission T (%)** : light intensity before and after cuvette

2.) **Absorbance**:  $Abs_{\lambda} = -\log(T)$  or „Extinction of light intensity“

3.) **Concentration**:

quantitative analysis of defined substance at a defined wavelength based on a calibration curve

# Transmission measurement



Transmission is the ratio of passed light  $I$  / initial light  $I_0$ :  $T \% = \frac{I}{I_0} \times 100$

Transmission measurement can be also used in photometers to measure **turbidity at 180°** (angle of light): => correction factor or turbidity measurement in the unit FAU



# Absorbance / Concentration Measurement

Absorbance = „Extinction of light“:  
Each substance has a specific spectra with an absorbance peak

=> Spectra is the preliminary task to define the wavelength for concentration measurement

Concentration measurement:

The specific wavelength at absorbance peak will be obtained by matching LED, optical filter or monochromator




# Relation %T / Abs / Conc.

## Transmission measurement:

The transmission of a sample varies **exponentially** with thickness and concentration

## Absorbance measurement:

Absorbance of a sample is **proportional** to thickness of the sample and concentration




Transmission (T%)	100	10	1	0,1	=> Logarithmic correlation
Absorbance $A = -\lg(T)$	0	1	2	3	} => linear correlation
Concentration (mg/l)	0	4	8	12	

# Concentration Measurement

Linearity for absorbance => concentration measurement based on

**calibration curve**

for each parameter (sustance):  
known concentrations at a defined  $\lambda$   
(at defined pathlength = cuvette size)



**Unknown concentration can be „read“ in the curve!**

Method data / Programs: programmed calibration curves in the meters  
=> automatic calculation of concentration by absorbance value



# Method data / Program for each parameter

Method data contain all necessary information for the concentration measurement:

$\lambda$  = usually absorbance peak of substance

Reagent blank  $E_0$  = coloration of reagent


Factor for calculation: slope of the curve

Citation & unit (e.g.  $\text{NO}_3\text{-N}$  mg/l)

Calculation factor to change citation & unit (e.g.  $\text{NO}_3$ ; mmol/l)


Sample blank (e.g. coloration of the sample) is **not** included!

„ $E_0$ “ would be increased: individually for sample, no defined factor!



# Pre-Requisites

- Coloured solution contains dissolved dye  
absorption of light leads to a coloration in **complementary color**  
Intensity is in relation to concentration
- The chemical reaction of parameter with reagent leads to **building or disappearance** (e.g. COD 4-40 mg/l) of dye in defined reaction time



Reaction must be selective for the substance – no reaction with other **disturbing substances**: see package leaflets

The developed dye must be stable for time of measurement => e.g. reading within 10 minutes after reaction time!

# WTW Photometer

## Portable Photometers

- pHotoFlex STD
- pHotoFlex pH
- pHotoFlex Turb




## Benchtop Photometers

- photoLab S6
- photoLab S12



## Spectrophotometer

- photoLab 7100
- photoLab 7600 UV-VIS



## **Poll Question #2**

What do you normally test using Lab and Field Photometers?

# pHotoFlex® Series Multiparameter Portable Colorimeter



# Water frame work directive – the portable Lab 2 Go!

## SERIES **pHoto Flex**

Environmental Protection:

### STD


Monitoring  
Comprehensive Water Analysis



Fish Farming and Surface Water:

### pH

Free Ammonia and CO<sub>2</sub>



Water, Beverage- and Quality Control:

### TURB

The real "Multi" with pH and Turbidity

### Accessories ...

... for Lab and on the Go



**LabStation**  
for even more comfort



**LSdata**-Software for

- GLP-compliant data management
- Programming of user-defined programs



**AQA:**  
pHotoFlex® Check



**Field Sets:**  
the mobile lab

**From Ammonia to Zinc:**

- More than 150 programs for routine analysis
- Test kits for every demand: from high precisely to easy, fast and economic
- Storage with Identification Number (ID)
- Datafilter for easy selection of data sets

**pHotoFlex® STD, plus:**

- Electrochemical pH-measurement from pH 0-16.00 (± 0.01) for Standard DIN pH-electrodes
- Automatic Temperature Compensation (ATC)
- Redox Measurement
- 1-3 point calibration with calibration interval setting and calibration protocol

**pHotoFlex® pH, plus:**

- IR light source acc. to DIN 27027/ISO 7027
- Turbidity measurement range from 0.01-1100 NTU/FNU providing drinking water accuracy
- Longterm stable AMCO Clear Standards with ±1% manufacturing accuracy
- The most versatile colorimeter for field and lab

# What PhotoFlex is commonly used for

## Disinfectants (e.g. Chlorine, Ozone)

- Fast degradation




## Environmental Test

- Nutrients (e.g.  $\text{NH}_4$ ,  $\text{NO}_3$ ,  $\text{PO}_4$ )



# photoLab S6 / S12

Routine measurement at it's best





# photoLab® S6/S12 – Just Measure!

Open lid, place cell, read result!



# photoLab® - S6/S12

## The best choice for routine tasks!

- Fastest measurement < 2 s to result
- Highest precision  
due to filter technology with ref.beam
- Most robust - no disalignment possible  
(„mechanics free“)

### photoLab Series offers:

- barcode for **round cells and reagent** test  
in rectangular cells 10, 20, 50 mm
- Automatic recognition for all cells



# What is S6/S12 commonly used for

- **Routine tests**
  - Dosing Applications
  - Regulatory Compliance
- **Repetitive tests**
- **Large Number of tests**  
e.g. Chemical Oxygen Demand,  
Nutrients, e.t.c.



# photoLab<sup>®</sup> 7000 Series



# photoLab<sup>®</sup> 7000 Series

## Routine and Spectral Analysis

### All in One – One for All

- Routine analysis with programmed method data
- Special and multistep procedures
- Comprehensive programming options
- Spectral and Kinetics Analysis



**From Water to Wine**  
**From R&D to Production**

# Extended PC Software

photoLab® color + photoLab Data *spectral*

- ✓ Color measurement from CIE to Gardner
- ✓ Data management offering link to LIMS (GLP Complaint)

Application: mostly Quality control

**Markets:**

Food&Beverage: Juice, Sugar, Oil  
Resin, varnish...




photoLab® color:  
Color measurement instead of color perception

- CIE color measurement: XYZ, x,y,z, CIE-L\* a\* b\*, CIE-L\* u\* v\*,
- Hazen (Pt-Co)
- Yellowness-Index
- ADMI color number
- ASTM
- Gardner
- Sugar color ICUMSA
- Beer color according to EBC and ASBC
- Iodine color number



# What is Photolab 7000 commonly used for

- **Where PhotoFlex or S6/S12 cannot fulfill.**
- **3<sup>rd</sup> Party Testing Lab**
  - E.g. Meet ISO 17025 standards
- **Teaching/Research**
  - E.g. Kinetics
- **Special Applications**
  - E.g. Sugars, Chlorophyll, Special color scales, e.t.c



# Takeaways

Think about how do you run your photometric tests?

## ❖ Field – Portability

- ✓ PhotoFlex

## ❖ Lab Routine – Ease of measurement from test to test

- ✓ Photolab S6/S12

➤ **Automated, stay tuned for next part of this presentation**

## ❖ Lab Special analysis/Research – Comprehensive functions

- ✓ Photolab 7000



# Instruments for Chemical Analysis




# Benjamin Chiang

- Lab and Process Business Development Manager
- 15 Years in process filtration and process instrumentation business
- 4 years with Xylem
- SEA Product manager for OIA analytical
- BDM for YSI, ISS & Aanderaa.  
MY,BN,PH



# Agenda

1. OIA Introduction & Portfolio
2. What is Automated Chemistry analyzer ( ACA)
3. Different ACA techniques and methods.
4. Advantages of automation
5. The FS3700
6. Application examples
7. Summary and closing



### **Poll Question #3**

What is your greatest concern when performing a sample analysis?



## Poll Question #4

How many samples of the same test (Eg. TKN, Cyanide) do you need to do per day?

# OI Analytical Laboratory and Online Products

## Total Organic Carbon (TOC) Analyzers ( Lab & Online)

- Detects and measures the concentration of organic matter and compounds in water samples.



## Gas Chromatography (GC) Products

- Selective GC detectors (PID, ELCD, FID, XSD, & PFPD) for GC instrumentation manufacturers, Purge & Trap (P&T) sample concentrators,

## Automated Chemistry Analyzers

- Continuous laboratory flow analyzers for water, soil, and plant samples.



## Cyanide Analyzer

- Measures available cyanide in precious metal leaching solutions by U.S. EPA Method





# FS3700 - Automating chemical analysis with OI Analytical



# Manual Method Process




Detector







# OIA Automated Process

Up to 360  
Samples in 1 run



Up to 2 types of chemical  
analysis in 1 run



Reagent

# Who Should Automate?

Anyone looking to improve laboratory efficiency

- Environmental
- Municipal Waste / Drinking water
- Agricultural
- Food and beverage
- Pharmaceutical
- Mining
- Petrochemical / Chemical
- Academic research
- Etc.



# Advantages of Automating your testings

Increases productivity.

Decreases labor cost.

Increases day to day reproducibility.

Reduces reagent use

Reduces the cost for waste disposal

Analysis costs are reduced.

All other relevant cost will be reduced



# Laboratory costs that can be reduced by automation

	<b>Laboratory Cost reduction by...</b>
<b>Direct labor</b>	<b>20 – 28 %</b>
<b>Indirect Labor</b>	<b>7 – 12 %</b>
<b>Operational Supplies</b>	<b>10 – 20 %</b>

# Two important automated chemistry technologies

## Automated Wet Chemistry

Continuous Flow

SFA


FIA

Robotic


Discrete



# Reagent usage is decreased by automation



# Waste generation is decreased by automation



# FS 3700 Automated Chemistry Analyzer Methods

Analyte	Technique	Method	Operating Range	MDL <sup>1</sup>	Throughput	Channel <sup>2</sup> Part #	Cartridge Part #
<b>Ammonia</b>	SFA, Gas Diffusion	USEPA 350.1	0.01-20.0 ppm 10.0-20,000 ppb	0.001 ppm 1.0 ppb	40 samples per hour	330109	330094
<b>Ammonia, Nitrogen</b> (Phenate)	FIA	USEPA 350.1	0.01-20 ppm	0.002 ppm	51 samples per hour	330353	330354
<b>Chloride</b>	SFA	Standard Methods 4500-Cl-E	1.0-200 ppm	0.12 ppm	60 samples per hour	330360	330361
<b>Cyanide</b> Available (1677)	FIA	OIA-1677-09	0.002-5.00 ppm 2.0-5,000 ppb	0.0005 ppm 0.5 ppb	30 samples per hour	330107	330092
<b>Cyanide</b> Available (D6888) (Sulfide abatement)	FIA	ASTM D6888-09	0.005-0.5 ppm 5.0-500 ppb	0.002 ppm 2.0 ppb	30 samples per hour	330106	330091
<b>Cyanide</b> Free (D7237)	FIA	ASTM D7237-10	2.0-500 ppb	0.5 ppb	30 samples per hour	330355	330356
<b>Cyanide</b> Free	Photometric Detection	ISO 14403	2.0-500 ppb	0.4 ppb	30 samples per hour	330371	330372
<b>Cyanide</b> Post-Distillation	FIA, Photometric Detection	USEPA 335.4	5.0-500 ppb	0.5 ppb	30 samples per hour	330351	330352
<b>Cyanide</b> Total	SFA, UV Digestion	ASTM D7511-09	0.003-0.5 ppm 3.0-500 ppb	0.0001 ppm 1.0 ppb	30 samples per hour	330076	330090
<b>Cyanide</b> Total	Photometric Detection	ISO 14403	2.0-500 ppb	0.4 ppb	30 samples per hour	330366	330367
<b>Hexavalent Chromium</b>	FIA	USEPA 600/4-79- 020	0.01-10 mg/L	0.0011 mg/L	48 samples per hour	331543	331544



Analyte	Technique	Method	Operating Range	MDL <sup>1</sup>	Throughput	Channel <sup>2</sup> Part #	Cartridge Part #
<b>MBAS</b>	Continuous Flow	ISO 16265	0.025-2.0 mg/L as LAS	0.008 mg/L as LAS	24 samples per hour	330357	330358
<b>Nitrate/Nitrite</b>	FIA	USEPA 353.2	0.01-10.0 ppm 10.0-10,000 ppb	0.001 ppm 1.0 ppb	60 samples per hour	330108	330093
	SFA		0.005-10.0 ppm	0.001 ppm	40 samples per hour	331377	331376
<b>Nitrate/Nitrite in Milk</b>	FIA w/ In-line Dialysis	ISO 14673-3	<b>Nitrate</b> 0.5 mg/L - 5.0 mg/L <b>Nitrite</b> 0.025 µg/L - 0.400 µg/L	<b>Nitrate</b> 0.016 mg/L <b>Nitrite</b> 0.0016 mg/L	30 samples per hour	331534	331535
<b>Phenol</b> In-line distillation	SFA	USEPA 420.2	1.0-500 ppb	0.5 ppb	22 samples per hour	330363	330364
<b>Phenol</b> Post-Distillation	FIA	USEPA 420.4	0.01-2.00 ppm 10.0-2,000 ppb	0.002 ppm 2.0 ppb	90 samples per hour	330110	330083
<b>Phosphorus</b> All Forms	FIA	USEPA 365.1	0.01-5.0 ppm 10.0-5,000 ppb	0.001 ppm 1.0 ppb	60 samples per hour	330111	330096
<b>Phosphorus</b> All Forms - Low Level	FIA	USEPA 365.1	0.001-0.1 ppm 1.0-100 ppb	0.0003 ppm 0.3 ppb	45 samples per hour	330112	330095
<b>Sulfate</b>	FIA Photometric	USEPA 375.2	1.0 mg/L - 25 mg/L	0.1 mg/L	40 samples per hour	331385	331386
<b>TKN</b> Total Kjeldahl Nitrogen	SFA, Gas Diffusion	USEPA 351.2	0.01-20.0 ppm 10.0-20,000 ppb	0.001 ppm 1.0 ppb	40 samples per hour	330109	330094

<sup>1</sup> Method Detection Limit (MDL) determined in accordance with 40 CFR Part 136 Appendix B

<sup>2</sup> Channels include the cartridge, detector, and valve (if required).

# Automated Chemistry Sample Types

We have Large Applications Library in these sectors. **ASK US!**

- Beverages
- Boiler feedwater
- Brackish waters
- Drinking water
- Feeds
- Fertilizers
- Foodstuffs
- Groundwater
- Wastewaters
- Industrial waste
- Mining process and wastewaters
- Plating baths
- Seawater
- Soil and plant extracts and digestion
- Surface waters
- Tobacco
- Wine
- Ultrapure water

(CFA) Continuous flow analyzers move liquid continuously through tubing.


**Within the tubes**

**Flow**

**Mix**


**React**

**Detect**






# (SFA) Segmented Flow Analysis minimizes dispersion with air bubbles.




(FIA) Flow injection analysis injects into an unsegmented stream.



# First Look – FS3700



Delivering VOC requested features




- Linear Logical Layout – Flows from left to right
- Modular hardware – add and modify on the fly
- Up to date software – no software issues
- Rapid and fast deployment – Easy to setup

# Chemistries available ( Top sellers )

- Available cyanide (EPA OIA-1677)
- Available cyanide (ASTM D6888)
- Total cyanide (ASTM D7511)
- Ammonia/TKN by gas diffusion
- Nitrate/Nitrite
- Phosphate, all forms
- Phosphate, all forms (low level)
- Phenol, Post-distillation



# What's in a channel?



- Everything needed to perform analysis is included in a channel.
- Cartridges include gas diffusion manifolds, heaters or UV lamps (as needed) all tubing, pump tubing and the chemistry kit.




# Everything is included?

Yes. The nesting structure of FS 3700 channels.



CHANNEL – Total Cyanide, ASTM D7511 330076



UV Ballast 330055 & Detector Cover 329515 are included at channel level.

# FlowView – user friendly interface



O.I. Analytical  
FlowSolution 3700 Analyzer Launcher

+

Add

×

Delete

⚙️

Modify

?

Help

▶

Training

10011010  
1101100  
011110

Upgrade

DEMO FS3700

O.I. Analytical  
Flow Solution 3700 Analyzer

Configure Analyzer

Configure System - Detectors

Primary Chassis

Sampler

Photometric CH1

Photometric CH2

Enable Detector 1

Mode: Photometric

Sample Gain: Photometric

Polarity: Amperometric

Cell Potential while in Standby mode (volts)

Enable Detector 2

Mode: Photometric

Sample Gain: 8 Reference Gain: 8

Polarity: Positive

Cell Potential while in Standby mode (volts)

Save

FS3700 State Machine: Idle Signal Poller: Idle Vial#: 0 Reqs: 0/0 State: None Active Sample Table: colortest

O.I. Analytical  
Flow Solution 3700 Analyzer

View Results

Channel 1

Channel 2

Channel 3

Channel 4

Summary

Print Results

Export Results

Save Results As

Open Results

30000

25000

20000

15000

10000

5000

0.0

0.0 200.0 400.0 600.0 800.0 1000.0 1200.0 1400.0 1600.0 1800.0 2000.0

Time (seconds)

Peak#	Cup	Sample Name	Rep	Type	Dil	Wt	Peak Height	Calc (ppb)	Avg Ht	Avg Calc (ppb)	%RSD	Flags
1	101	SYNC_100 ppb	1	SPL	1		29117.0	0.0	29117.0	0.0	0.0	
2	102	Blank	1	SPL	1		33.0	0.0				
3	102	Blank	2	SPL	1		4.0	0.0	18.5	0.0	0.0	
4	103	0 ppm PO4	1	STD1	1		4.0	0.0	4.0	0.0	0.0	
5	104	1 ppb PO4	1	STD2	1		226.0	1.0	226.0	1.0	0.0	
6	105	5 ppb PO4	1	STD3	1		1371.0	5.0	1371.0	5.0	0.0	
7	106	10 ppb PO4	1	STD4	1		2885.0	10.0	2885.0	10.0	0.0	
8	107	25 ppb PO4	1	STD5	1		7318.0	25.0	7318.0	25.0	0.0	
9	108	50 ppb PO4	1	STD6	1		14683.0	50.0	14683.0	50.0	0.0	

Result File: 2014-06-05\_LL\_PO4\_Linearity\_Accuracy

# Flow Solution 3700

Some of the applications recently used...

- Monterrey Bay Analytical Services. Environmental application. FS3700 Purchase with Ammonia/TKN and Phosphorus (All Forms).
- Gerhardt Block Digestor (for TKN or TP digestions)
- KJELDATHERM Digest System KBL40S 230V coupled with the FS3700



KJELDATHERM®



# Flow Solution S3700

## Another application recently used...

- Unique application for the FS3700 is Nitrite/Nitrate in Milk Products
- ISO Method 14673 (Flow Injection) ISO Method 14673-2 (Segmented)
- Determination of NO<sub>2</sub>/NO<sub>3</sub> contents of milk and milk products by FIA. Applicable to hard, semi-hard, soft cheeses, whey powder, milk powder and milk based infant foods.
- Tillamook in Boardman, OR (Cheese)
- Dairy America in Tulare, CA



# Flow Solution S3700

Another application recently used...

- Unique application for the FS3700 is Cyanide Analysis in Gold Mining
- USEPA Methods OIA-1677, OIA-1678 & Total CN Method ASTM D7511-09
- Determination of cyanide levels in tailings ponds of Gold mines
- Determination of cyanide levels released into the environment
- Mines in USA, Africa(Ghana) and Australia.



# Summary on Flow Solution 3700

- The FS 3700 introduces **flexible and innovative technology** to automated your chemistry smoothly.
- Icon-driven software that is as **easy to use** like a smartphone
- **Modular**, plug-n-play components allows **easy customization** on the system for your chemistries and improve your laboratory workflow.
- **Simple** to setup, operate and **inexpensive** to maintain.
- **High sample throughput**, while handling complex matrices that other analyzers can't process.
- **Minimal Operator involvement**. Once set, move to other jobs...
- FS3700 handles a **wide variety of chemistries** such as Total and Free Cyanide, NO<sub>2</sub>/NO<sub>3</sub>, Ammonia/TKN, TP and low-level Ortho-P.



## **Poll Question #5**

Would you like someone from Xylem to contact you about photometry or automated chemical analyzer solutions?

# Questions?

**Contact us:**

**Xavier Tan (Photometry)**

[Xavier.Tan@xyleminc.com](mailto:Xavier.Tan@xyleminc.com)

**Benjamin Chiang (Auto Chemical Analyzers)**

[Benjamin.Chiang@xyleminc.com](mailto:Benjamin.Chiang@xyleminc.com)